




# Intubation Aids



## Tube Guides

### ENDOTRACHEAL TUBE INTRODUCERS (BOUGIE)

- Sterile, disposable & latex free
- Manufactured from low density polyethylene (provides proper stiffness for ease of insertion)
- Calibrated (distance of insertion easily observed for safety)
- Coude tip (facilitates insertion in adults)
- Pediatric size fits in 4mm to 6mm tubes
- Adult size fit in 6mm to 11mm tubes
- Pedi sizes have calibrated markings on both ends allowing it to be used as a straight or coude tip




CAT#	DESCRIPTION	BOX
9-0211-70	Pedi 10FR x 70cm w/Coude tip	10
9-0211-71	Pedi 10FR x 70cm w/Coude tip	100
9-0212-70	Adult 15FR x 70cm w/Coude tip	10
9-0212-71	Adult 15FR x 70cm w/Coude tip	100
9-0212-72	Adult 15FR x 70cm w/Straight tip	10
9-0212-73	Adult 15FR x 70cm w/Straight tip	100



### MALLEABLE INTRODUCING STYLETS

Malleable introducing stylet contain an internal wire which provides the user the ability to pre-bend the stylet, thus controlling the direction when inserting endotracheal tubes, bronchial catheters and laryngeal type masks.

CAT#	DESCRIPTION	BOX
9-0213-82	Pedi Malleable 10FR x 70cm w/Coude & Straight tip	10
9-0213-92	Adult Malleable 15FR x 70cm w/Coude tip	10



### DISPOSABLE ALUMINUM STYLETS

- Malleable with rounded end
- Sterile & latex free
- Satin finish to reduce friction
- Individually packaged
- Pediatric stylet fits in 2-6mm tubes
- Adult stylet fits in 6-11mm tubes

CAT#	DESCRIPTION	BOX
9-0204-11	Disp. Pediatric Aluminum Stylet, 6FR, 11"	25
9-0204-25	Disp. Adult Aluminum Stylet, 9FR, 16"	25



# Intubation Aids



## Stylets

CAT#	DESCRIPTION	BOX
9-0204-00	Adult Aluminum Stylet, 10FR, 16"	1
9-0204-10	Adult Aluminum Stylet, 10FR, 16"	10
9-0204-12	Pedi Aluminum Stylet, 10FR, 12"	1
9-0204-22	Pedi Aluminum Stylet, 10FR, 12"	10

**Single Patient Use**



CAT#	DESCRIPTION	BOX
9-0207-03	Surgical Stainless Steel, <b>Reusable</b> Adult, 6FR, 16.5"	1
9-0207-01	Surgical Stainless Steel, <b>Reusable</b> Child, 6FR, 11.5"	1



CAT#	DESCRIPTION	BOX
9-0206-01	Hollow Copper Tube with Suction Port 16.5"	1

**Single Patient Use**




CAT#	DESCRIPTION	BOX
9-0203-00	Solid Copper Malleable with Adjustable Stop, 8FR, 14.75"	1

**Single Patient Use**



CAT#	DESCRIPTION	BOX
9-0204-18	Adult Fiber Optic Stylet with Adjustable Stop 18" <b>Reusable</b>	1
9-0204-14	Pedi Fiber Optic Stylet with Adjustable Stop 14" <b>Reusable</b>	1
9-0204-19	Adult Replacement Bundle	1
9-0204-15	Pedi Replacement Bundle	1
5-0240-59	Lamp for Fiber Optic Stylet	1
9-0210-05	Replacement Stops	5





# Intubation Aids



## Surgical Instruments

### MAGILL FORCEPS CLOSED-TIP

- Surgical stainless steel

CAT#	DESCRIPTION	BOX
9-3018-01	Infant Magill Closed-Tip Forceps	1
9-3018-02	Child Magill Closed-Tip Forceps	1
9-3018-03	Adult Magill Closed-Tip Forceps	1



### MAGILL FORCEPS OPEN-TIP

- Surgical stainless steel

CAT#	DESCRIPTION	BOX
9-3016-02	Child Magill Open-Tip Forceps	1
9-3016-03	Adult Magill Open-Tip Forceps	1



### AILLON TUBE BENDER

One of the most common problems associated with nasotracheal intubation is observed at a time when the endotracheal tube is about to cross the vocal cords but stops making progress. The tube tends to follow a course which is forward and upward. The trachea's main direction is downward. The solution is to have a tool that not only will control the tip of the tube, but will also bend it down in order to conform with the direction of the trachea. The Aillon forceps tube bender is the solution.

CAT#	DESCRIPTION	BOX
9-3021-00	Aillon Tube Bender	1



### TONGUE SEIZING FORCEPS

Made of surgical stainless steel with latex rubber tips.

CAT#	DESCRIPTION	BOX
9-3021-50	Tongue Seizing Forceps	1

